## **Safety Needles for Sharp**

NEEDLE	TIP / DIAMETER	CHARACTERISTICS	PMU TREATMENT	PMU TECHNIQUE
1-point nano	Round 0,25 mm	Our thinnest needle creates ultra-sharp, crispy lines	Eyebrows, full lips, lip liner and areola capillary	Hairstroke, lining, pointilisme, shading, one way striping, two way striping
1-point liner	Round 0,40 mm	Suitable for detailed lines	Eyebrow shading, eyeliner and lip liner	One way striping, pointilism, hairstroke, shading
1-point nano tube	Round 0,25 mm	Ultra-fine, high quality needle	Eyebrows, full lips, lip liner and areola capillary	Hairstroke, lining, one way striping, two way striping, pointilisme, shading
3-point micro	Round 0,18 mm	Suitable for ultra-fine lines	Eyebrows, eyeliner, eyelash enhancement, eyeshadow, full lips, glossy lips, lip liner, areola and scars	One way striping, two way striping, lining, shading
3-point outline	Round 0,25 mm	One of the points protrudes in the direction of the supply opening	Eyebrows, eyeliner, eyelash enhancement, full lips and lip liner	One way striping, pointilism, lining, shading
3-point nano sloped	Round 0,25 mm	Flat form comprising 3 ultra-thin needles in a row	Eyebrows, eyeliner, eyelash enhancement, eyeshadow, lip contour, full lips and areola	Hairstroke, shading, lining
5-point magnum	Round 0,30 mm	Flat form comprising 5 needles in 2 rows	Eyeliner, full lips, areola and scars	One way striping, two way striping, circling, lining, shading